

Product Datasheet

Human IL1RAP ELISA Kit (orb1201935)

Description	Human Interleukin 1 Receptor Accessory Protein (IL1RAP) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q9NPH3
Sensitivity	0.065ng/mL
Entrez	3556
Expiration Date	Please enquire.