

Product Datasheet

Human FGF12 ELISA Kit (orb1201958)

Description	Human Fibroblast Growth Factor 12 (FGF12) ELISA Kit
Reactivity	Human
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P61328
Sensitivity	8.1pg/mL
Entrez	2257
Expiration Date	Please enquire.