

Product Datasheet

Human FGF11 ELISA Kit (orb1201959)

Description	Human Fibroblast Growth Factor 11 (FGF11) ELISA Kit
Reactivity	Human
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q92914
Sensitivity	5.9pg/mL
Entrez	2256
Expiration Date	Please enquire.