

Product Datasheet

Human C1QTNF3 ELISA Kit (orb1201961)

Description	Human C1q And Tumor Necrosis Factor Related Protein 3 (C1QTNF3) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q9BXJ4
Sensitivity	0.127ng/mL
Entrez	114899
Expiration Date	Please enquire.