

Product Datasheet

Human ATP4a ELISA Kit (orb1201965)

Description	Human ATPase, H ⁺ /K ⁺ Exchanging Alpha Polypeptide (ATP4a) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P20648
Sensitivity	0.055ng/mL
Entrez	495
Expiration Date	Please enquire.