

Product Datasheet

Rat TGFb1 ELISA Kit (orb1201983)

Description	Rat Transforming Growth Factor Beta 1 (TGFb1) ELISA Kit
Reactivity	Rat
Range	31.2-2000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, platelet-poor plasma, tissue homogenates, cell culture supernates or other biological fluids
Uniprot ID	P17246
Sensitivity	13.3pg/mL
Entrez	59086
Expiration Date	Please enquire.