

Product Datasheet

Rat TGFa ELISA Kit (orb1201985)

Description	Rat Transforming Growth Factor Alpha (TGFa) ELISA Kit
Reactivity	Rat
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P01134
Sensitivity	6.7pg/mL
Entrez	24827
Expiration Date	Please enquire.