

Product Datasheet

Chicken IL4 ELISA Kit (orb1202011)

Description	Chicken Interleukin 4 (IL4) ELISA Kit
Reactivity	Gallus
Range	7.81-500pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q5W4U1
Sensitivity	2.7pg/mL
Entrez	416330
Expiration Date	Please enquire.