

Product Datasheet

Chicken IL1b ELISA Kit (orb1202025)

Description	Chicken Interleukin 1 Beta (IL1b) ELISA Kit
Reactivity	Gallus
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	073909
Sensitivity	5.7pg/mL
Entrez	395196
Expiration Date	Please enquire.