

Product Datasheet

Rabbit IL1a ELISA Kit (orb1202026)

Description	Rabbit Interleukin 1 Alpha (IL1a) ELISA Kit
Reactivity	Rabbit
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P04822
Sensitivity	5.4pg/mL
Entrez	100009250
Expiration Date	Please enquire.