

Product Datasheet

Rat IL15 ELISA Kit (orb1202036)

Description	Rat Interleukin 15 (IL15) ELISA Kit
Reactivity	Rat
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P97604
Sensitivity	5.4pg/mL
Entrez	25670
Expiration Date	Please enquire.