

Product Datasheet

Rat IL11 ELISA Kit (orb1202043)

Description	Rat Interleukin 11 (IL11) ELISA Kit
Reactivity	Rat
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q99MF5
Sensitivity	5.8pg/mL
Entrez	171040
Expiration Date	Please enquire.