

Product Datasheet

Mouse IL10Ra ELISA Kit (orb1202045)

Description	Mouse Interleukin 10 Receptor Alpha (IL10Ra) ELISA Kit
Reactivity	Mouse
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q61727
Sensitivity	29pg/mL
Entrez	16154
Expiration Date	Please enquire.