

Product Datasheet

Mouse IGFBP7 ELISA Kit (orb1202048)

Description	Mouse Insulin Like Growth Factor Binding Protein 7 (IGFBP7) ELISA Kit
Reactivity	Mouse
Range	6.25-400pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q61581
Sensitivity	2.34pg/mL
Entrez	29817
Expiration Date	Please enquire.