

Product Datasheet

Rat IGFBP1 ELISA Kit (orb1202054)

Description	Rat Insulin Like Growth Factor Binding Protein 1 (IGFBP1) ELISA Kit
Reactivity	Rat
Range	1.56-100ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P21743
Sensitivity	0.55ng/mL
Entrez	25685
Expiration Date	Please enquire.