

Product Datasheet

Chicken IFN γ ELISA Kit (orb1202060)

Description	Chicken Interferon Gamma (IFN γ) ELISA Kit
Reactivity	Gallus
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P49708
Sensitivity	5.8pg/mL
Entrez	396054
Expiration Date	Please enquire.