

Product Datasheet

Rat EG-VEGF ELISA Kit (orb1202081)

Description	Rat Endocrine Gland Derived Vascular Endothelial Growth Factor (EG-VEGF) ELISA Kit
Reactivity	Rat
Range	31.2-2000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q8R414
Sensitivity	13.6pg/mL
Entrez	192205
Expiration Date	Please enquire.