

Product Datasheet

Rabbit EDN1 ELISA Kit (orb1202084)

Description	Rabbit Endothelin 1 (EDN1) ELISA Kit
Reactivity	Rabbit
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	competitive inhibition
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P29560
Sensitivity	5.63pg/mL
Entrez	100009270
Expiration Date	Please enquire.