

Product Datasheet

Rat BMP4 ELISA Kit (orb1202090)

Description	Rat Bone Morphogenetic Protein 4 (BMP4) ELISA Kit
Reactivity	Rat
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q06826
Sensitivity	6.1pg/mL
Entrez	25296
Expiration Date	Please enquire.