

Product Datasheet

Mouse HSPb6 ELISA Kit (orb1202117)

Description	Mouse Heat Shock Protein Beta 6 (HSPb6) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q5EBG6
Sensitivity	0.067ng/mL
Entrez	243912
Expiration Date	Please enquire.