

Product Datasheet

Mouse CASP8 ELISA Kit (orb1202124)

Description	Mouse Caspase 8 (CASP8) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	089110
Sensitivity	0.069ng/mL
Entrez	12370
Expiration Date	Please enquire.