

Product Datasheet

Human IL2Rg ELISA Kit (orb1202164)

Description	Human Interleukin 2 Receptor Gamma (IL2Rg) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P31785
Sensitivity	0.067ng/mL
Entrez	3561
Expiration Date	Please enquire.