

Product Datasheet

Human IL26 ELISA Kit (orb1202167)

Description	Human Interleukin 26 (IL26) ELISA Kit
Reactivity	Human
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9NPH9
Sensitivity	6.1pg/mL
Entrez	55801
Expiration Date	Please enquire.