

Product Datasheet

Human IL13Ra2 ELISA Kit (orb1202176)

Description	Human Interleukin 13 Receptor Alpha 2 (IL13Ra2) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q14627
Sensitivity	0.129ng/mL
Entrez	3598
Expiration Date	Please enquire.