

Product Datasheet

Human FGFR1 ELISA Kit (orb1202188)

Description	Human Fibroblast Growth Factor Receptor 1 (FGFR1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P11362
Sensitivity	0.057ng/mL
Entrez	2260
Expiration Date	Please enquire.