

Product Datasheet

Human FGF10 ELISA Kit (orb1202189)

Description	Human Fibroblast Growth Factor 10 (FGF10) ELISA Kit
Reactivity	Human
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	015520
Sensitivity	5.7pg/mL
Entrez	2255
Expiration Date	Please enquire.