

Product Datasheet

Human CSF2Rb ELISA Kit (orb1202192)

Description	Human Colony Stimulating Factor 2 Receptor Beta (CSF2Rb) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P32927
Sensitivity	0.061ng/mL
Entrez	1439
Expiration Date	Please enquire.