

Product Datasheet

Human FABP2 ELISA Kit (orb1202223)

Description	Human Fatty Acid Binding Protein 2, Intestinal (FABP2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P12104
Sensitivity	0.057ng/mL
Entrez	2169
Expiration Date	Please enquire.