

## Product Datasheet

### Human CA19-9 ELISA Kit (orb1202234)

|                            |                                                                                                      |
|----------------------------|------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Human Carbohydrate Antigen 19-9 (CA19-9) ELISA Kit                                                   |
| <b>Reactivity</b>          | Human                                                                                                |
| <b>Range</b>               | 3.12-200U/mL                                                                                         |
| <b>Tested Applications</b> | ELISA                                                                                                |
| <b>Note</b>                | For research use only                                                                                |
| <b>Application notes</b>   | sandwich                                                                                             |
| <b>Sample Types</b>        | serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| <b>Sensitivity</b>         | 1.18U/mL                                                                                             |
| <b>Expiration Date</b>     | Please enquire.                                                                                      |