

Product Datasheet

Mouse TGFa ELISA Kit (orb1202254)

Description	Mouse Transforming Growth Factor Alpha (TGFa) ELISA Kit
Reactivity	Mouse
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P48030
Sensitivity	5.5pg/mL
Entrez	21802
Expiration Date	Please enquire.