

Product Datasheet

Mouse IGFBP4 ELISA Kit (orb1202275)

Description	Mouse Insulin Like Growth Factor Binding Protein 4 (IGFBP4) ELISA Kit
Reactivity	Mouse
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P47879
Sensitivity	6.4pg/mL
Entrez	16010
Expiration Date	Please enquire.