

Product Datasheet

Mouse FGF4 ELISA Kit (orb1202286)

Description	Mouse Fibroblast Growth Factor 4 (FGF4) ELISA Kit
Reactivity	Mouse
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P11403
Sensitivity	5.4pg/mL
Entrez	14175
Expiration Date	Please enquire.