

Product Datasheet

Mouse CTGF ELISA Kit (orb1202293)

Description	Mouse Connective Tissue Growth Factor (CTGF) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P29268
Sensitivity	0.061ng/mL
Entrez	14219
Expiration Date	Please enquire.