

Product Datasheet

Human IFNa2 ELISA Kit (orb1202346)

Description	Human Interferon Alpha 2 (IFNa2) ELISA Kit
Reactivity	Human
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P01563
Sensitivity	6.3pg/mL
Entrez	3440
Expiration Date	Please enquire.