

Product Datasheet

Human FIL1d ELISA Kit (orb1202348)

Description	Human Interleukin 1 Delta (FIL1d) ELISA Kit
Reactivity	Human
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9UBH0
Sensitivity	5.7pg/mL
Entrez	26525
Expiration Date	Please enquire.