

Product Datasheet

Human FGF6 ELISA Kit (orb1202350)

Description	Human Fibroblast Growth Factor 6 (FGF6) ELISA Kit
Reactivity	Human
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P10767
Sensitivity	37.33pg/mL
Entrez	2251
Expiration Date	Please enquire.