

Product Datasheet

Human ADAM17 ELISA Kit (orb1202367)

Description	Human A Disintegrin And Metalloprotease 17 (ADAM17) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P78536
Sensitivity	0.056ng/mL
Entrez	6868
Expiration Date	Please enquire.