

Product Datasheet

Rat FGF2 ELISA Kit (orb1202377)

Description	Rat Fibroblast Growth Factor 2, Basic (FGF2) ELISA Kit
Reactivity	Rat
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P13109
Sensitivity	4.37pg/mL
Entrez	54250
Expiration Date	Please enquire.