

Product Datasheet

Human SPD ELISA Kit (orb1202407)

Description	Human Surfactant Associated Protein D (SPD) ELISA Kit
Reactivity	Human
Range	6.25-400ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, lung lavage fluid, cell culture supernates or other biological fluids.
Uniprot ID	P35247
Sensitivity	2.16ng/mL
Entrez	6441
Expiration Date	Please enquire.