

Product Datasheet

Guinea TNFa ELISA Kit (orb1202420)

Description	Guinea pig Tumor Necrosis Factor Alpha (TNFa) ELISA Kit
Reactivity	Guinea pig
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P51435
Sensitivity	5.5pg/mL
Entrez	100135630
Expiration Date	Please enquire.