

Product Datasheet

Human KLK10 ELISA Kit (orb1202425)

Description	Human Kallikrein 10 (KLK10) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	O43240
Sensitivity	0.048ng/mL
Entrez	5655
Expiration Date	Please enquire.