

Product Datasheet

Guinea IL1b ELISA Kit (orb1202427)

Description	Guinea pig Interleukin 1 Beta (IL1b) ELISA Kit
Reactivity	Guinea pig
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9WVG1
Sensitivity	25pg/mL
Entrez	100135556
Expiration Date	Please enquire.