

Product Datasheet

Human MMP7 ELISA Kit (orb1202450)

Description	Human Matrix Metalloproteinase 7 (MMP7) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P09237
Sensitivity	0.055ng/mL
Entrez	4316
Expiration Date	Please enquire.