

Product Datasheet

Bovine FN ELISA Kit (orb1202452)

Description	Bovine Fibronectin (FN) ELISA Kit
Reactivity	Bovine
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P07589
Sensitivity	0.066ng/mL
Entrez	280794
Expiration Date	Please enquire.