

## Product Datasheet

### Rabbit TNFa ELISA Kit (orb1202505)

|                            |                                                                                                      |
|----------------------------|------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit Tumor Necrosis Factor Alpha (TNFa)<br>ELISA Kit                                               |
| <b>Reactivity</b>          | Rabbit                                                                                               |
| <b>Range</b>               | 15.6-1000pg/mL                                                                                       |
| <b>Tested Applications</b> | ELISA                                                                                                |
| <b>Note</b>                | For research use only                                                                                |
| <b>Application notes</b>   | sandwich                                                                                             |
| <b>Sample Types</b>        | serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| <b>Uniprot ID</b>          | <b>P04924</b>                                                                                        |
| <b>Sensitivity</b>         | 5.6pg/mL                                                                                             |
| <b>Entrez</b>              | <b>100009088</b>                                                                                     |
| <b>Expiration Date</b>     | Please enquire.                                                                                      |