

Product Datasheet

Rabbit IL1b ELISA Kit (orb1202508)

Description	Rabbit Interleukin 1 Beta (IL1b) ELISA Kit
Reactivity	Rabbit
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P14628
Sensitivity	33pg/mL
Entrez	100008990
Expiration Date	Please enquire.