

Product Datasheet

Human IL1z ELISA Kit (orb1202539)

Description	Human Interleukin 1 Zeta (IL1z) ELISA Kit
Reactivity	Human
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9NZH6
Sensitivity	24pg/mL
Entrez	27178
Expiration Date	Please enquire.