

Product Datasheet

Human VEGFD ELISA Kit (orb1202561)

Description	Human Vascular Endothelial Growth Factor D (VEGFD) ELISA Kit
Reactivity	Human
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	O43915
Sensitivity	6.1pg/mL
Entrez	2277
Expiration Date	Please enquire.