

Product Datasheet

Canine VEGFA ELISA Kit (orb1202567)

Description	Canine Vascular Endothelial Growth Factor A (VEGFA) ELISA Kit
Reactivity	Canine
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9MYV3
Sensitivity	5.4pg/mL
Entrez	403802
Expiration Date	Please enquire.