

Product Datasheet

Human ITaC ELISA Kit (orb1202580)

Description	Human Interferon Inducible T-Cell Alpha Chemoattractant (ITaC) ELISA Kit
Reactivity	Human
Range	31.2-2000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	O14625
Sensitivity	10.8pg/mL
Entrez	6373
Expiration Date	Please enquire.